



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **PHILIPS BRILLIANCE 64 computer tomograph (Reconditioned)**

Used, very good technical condition,  
Dutch production,  
Year of manufacture: 2006,  
64-row,  
Power: 60 kW,  
Frequency: 50/60 Hz,  
Segmented anode: the division of the anode into twelve segments compensates for the adverse effects of heating and cooling cycles, increasing system reliability,  
Dynamic focus movement doubles the number of projections and increases spatial resolution in the imaging plane,  
Gantry hole diameter: 700 mm,  
Intercom: two-way connection between the gantry and the console,  
Gantry tilt angle: -30 ° to + 30 °, in 0.5 ° increments,  
Focus to detector distance: 1040 mm,  
Focus-isocenter distance: 570 mm,  
Vertical movement range: from 578 to 1028 mm, adjustable every 1.0 mm,  
Manual table travel: 1900 mm,  
Scanning range: 1750 mm,  
Table load: 204 kg,  
Carbon fiber table top with foot pedal and manual controls for easy positioning and quick release  
Anode heat capacity: 8 MHU,  
Anode diameter: 200 mm,  
Tube type: PHILIPS MRC 800,  
Year of tube production: 2015,  
Number of scans: 854 846 scan sec.,  
Software version: 2.6.2.21004,  
Console type: Dell 670,  
Detector: (UDMS),  
Has a valid Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (only Poland): Installments, Leasing, Loan,

