



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table EOS M04 (New)

New,
Made in Poland,
Sturdy steel structure coated with powder paint,
Standard features include three-section sleep surface, foldable to Pivot and seat positions, made of plywood,
Height adjustment of the table is carried out with a hand crank at the back of the table frame,
Fully ball-bearing mounted height adjustment mechanism,
Individual parts of the tabletop are adjusted by means of gas springs,
There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available at an extra charge,
Table features and benefits
Gas spring from the German company BANSBACH,
Upholstery of very good quality,
Foam, 5 cm thick,
Steel ball bearings are used in the movable parts of the table, which makes the table very reliable,
Frame elements made of the best Polish steel,
As additional options are available (for an extra charge):
Armrests,
Side shelves,
Hanger for the underlay,
Castor system,
Euro rails (40cm),
Knee pads,
Injection support,
Arm support for the chair,
Technical parameters:
Length 200cm,
Width 65/70/80 cm,
Height 48 - 100 cm,
Headrest angle adjustment -700 to 400,
Backrest angle adjustment 00 to 850,
Adjustment of the angle of the top 00 to 300,
Own weight 90 kg,
Maximum load 160 kg,
It has a current Technical Passport,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE



