



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table EOS M05 (New)

New,
Made in Poland,
Stable steel construction coated with powder paint,
Standard features one-piece plywood sleep surface,
Height adjustment of the table is carried out with a hand crank placed in the back part of the table frame,
The height adjustment is carried out with a hand crank placed in the back part of the table frame, the mechanism of height adjustment is fully bearing-mounted,
A characteristic element of the table is a wide tabletop, which enables to conduct exercises with the use of the Vojta and Bobath methods,
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours, which are available for an extra charge,
Additional options are available (at an extra charge):
Sleeper hanger,
Mobile castor system,
Technical parameters:
Length 200 cm,
Width from 80 to 120 cm,
Height 48 - 100 cm,
Weight 120 kg,
Maximum load 160 kg,
It has a current Technical Passport,
Construction warranty 24 months,
Warranty on hydraulic accessories 12 months,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

