



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



Rehabilitation table EOS H05 (New)

New,
 Made in Poland,
 It has a stable steel construction coated with powder paint,
 Standard features include one-piece plywood sleep surface,
 Height adjustment of the table is carried out with the use of a hydraulic actuator by repeatedly pressing with your foot the lever located on the actuator,
 The height adjustment is carried out with a hydraulic cylinder by repeatedly pressing the lever at the cylinder with your foot,
 A characteristic element of the table is a wide tabletop, which enables to conduct exercises with the use of the Vojta and Bobath methods,
 It is possible to order upholstery in the colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,
 Features and advantages of the table
 Upholstery of very good quality,
 The foam is 5 cm thick,
 Steel ball bearings are used in the movable parts of the table, which ensures unreliability of the table,
 Frame elements made of the best Polish steel,
 As additional options are available (for an extra charge):
 Hanger for the underlay,
 Mobile castor system,
 Technical parameters:
 Length 200 cm,
 Width from 80 to 120 cm,
 Height 48 - 100 cm,
 Self-supporting weight 124 kg,
 Maximum load 160 kg,
 It has a current Technical Passport,
 Warranty on structure 24 months,
 Warranty on hydraulic accessories 12 months,
 If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH
 HABLAMOS ESPANOL
 WIR SPRECHEN DEUTSCH
 ON PARLE FRANCAIS
 WIJ SPREKEN NEDERLANDS
 МЫ ГОВОРИМ ПО РУССКИ
 МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
 WE SHIP WORLDWIDE

