



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table ALLOS M (New)

New,
Made in Poland,
Stable steel construction powder coated,
Height adjustment of the table is carried out with a hand crank placed in the back part of the table frame,
Standard features include two-part plywood sleep surface,
Headrest is adjustable with a gas spring,
There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available at an extra charge,
Table features and benefits
Gas spring from the German company BANSBACH,
Upholstery of very good quality,
Foam, 5 cm thick,
Steel ball bearings are used in the movable parts of the table, which makes the table very reliable,
Frame elements made of the best Polish steel,
As additional options are available (for an extra charge):
Armrests,
Hanger for pad,
Castor system,
Euro rails (40cm),
Knee pads,
Injection support,
Technical parameters:
Length 200cm
Width 65/70/80 cm
Height 52-90 cm
Headrest angle adjustment -70° to 40°
82 kg empty weight
Maximum load 160 kg
It has a valid Technical Passport,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

