



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table ASES E (None)

New,
Made in Poland,
Sturdy steel structure coated with powder paint,
Standard features include a three-section plywood sleep surface,
The central part of the bed can be raised to the angle of 30°, which facilitates USG examinations,
Height adjustment of the table is carried out with an electric actuator, in this model a fully bearing-mounted height adjustment mechanism is used,
Individual parts of the tabletop are adjustable with gas springs,
There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available at an extra charge,
Table features and benefits
Gas spring from German company BANSBACH,
Electric actuator by a German company DEWERT,
Upholstery of very good quality,
5 cm thick foam,
Steel ball bearings were used in the movable parts of the table, which makes the table very reliable,
Frame elements made of the best Polish steel,
As additional options are available (for an extra charge):
Armrests,
Hanger for pad,
System of castors,
Foot remote control,
Control from the frame,
Euro rails (40cm),
Knee bolsters,
Injection support,
Technical parameters:
Length 200cm
Width 65/70/80 cm
Height 48-100 cm
Headrest angle adjustment 50° to - 70°
Adjustment of middle part angle 0° to 30°
own weight 92 kg
Maximum load 160 kg
It has a valid Technical Passport,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE



