



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Diagnostic and treatment table URBAN (Electric) (None)

New,

Polish production,

It has a stable steel structure covered with powder paint,

Professional table for Chiropractic, Osteopathy, Manual Therapy and Therapeutic Massage,

The model of the offered table has functions that are necessary to perform therapy in lean, obese as well as pregnant women,

The advantage of the table is the possibility of lifting with the help of gas springs, separately or together the lumbar part and the chest part in the case of obese people or pregnant women,

An important function in the offered model is the wide adjustment range of the headrest. A modern solution for adjusting the position up and down with the inclination angle setting, allows for comfortable placement of the patient with small and large chest and women with large breasts,

Adjusting the width of the headrest elements will help adjust the parameters to the shape of the face. A wide range of adjustments makes the patient not dependent on slim or obese will feel comfortable,

Strong table construction ensures maximum stability while working with the patient, and the applied solutions and functions as well as quiet and smooth work guarantee maximum comfort for the staff and the patient,

The height of the table is adjusted by means of an electric actuator controlled by a hand-held remote control,

There is a free opportunity to order upholstery in the color of your choice, except for the colors 'PRESTIGE' which are available at extra cost,

As additional options are available (at extra cost):

Hanger for paper,

Frame control,

Foot control,

A system of castors,

Technical data:

Length: 195 cm,

Width 65/70/80 cm,

Height: 48-100 cm,

Adjustment of the headrest: - 30 ° to 30 °,

Height adjustable headrest: 26 cm,

Slit headrest: From 6 to 12 cm,

Adjustment of the backrest angle: 0 ° to 15 °,

Adjustment of the angle of the central part: 0 ° to 25,

Own weight 93 kg,

Maximum load: 160 kg,

Electric power supply: Voltage: 230V / 50-60 Hz, power: max160W, protection class: II, housing: IP54,

Has a current Technical Passport issued,

Warranty for 24 months constructions,

Warranty for hydraulic and electric accessories 12 months

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



