



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion pump (syringe) BECTON DICKINSON PILOT C (Reconditioned)

Used, very good technical condition,
French production,
Power supply: 230 V and through a built-in 6 V / 1.1 Ah battery,
Frequency: 50/60 Hz,
Power: 23 VA,
Readable electronic graphic display and additional liquid crystal display of line pressure,
Syringe selection: 20 and 50/60 ml,
The rate of supply: 0.1 - 120 ml / h,
Bolus: 0.1 - 800 ml / h,
Infusion volume: 0.1 - 999.9 ml,
Acoustic and visual alarms,
Syringe size reader,
Programmable drug name,
Automatic bolus reduction after occlusion,
Infusion pressure display and settings,
Quick start,
RS232 socket,
Convenient carrying handle,
Dimensions: 32.5 x 19 x 12 cm,
Weight: 2.6 kg,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Installments, Leasing, Loan,

