



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STERIS QRC1752INT endoscopic connection (Reconditioned)

Used, very good technical condition,
U.S.A. production,

The QRC1752INT endoscopic terminal is intended only for use with the SYSTEM 1 endoscope washer and the C1140S1 container / tray and for the SYSTEM 1 EXPRESS or SYSTEM 1 PLUS washer with the C1140INT tray,

The connector is compatible with RICHARD WOLF 7326071 and 7326076 endoscopes,

The set includes a package of 100 disposable filters,

Multi-language user manual,

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month on the international market,

Financing option (Poland only): Installments, Leasing, Loans,

