



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical laser PHOTOMEDEX LASERPRO 980 25W (Reconditioned)

Used, perfect technical condition,
Production U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power consumption: 540 VA,
Laser type: semiconductor diode laser,
Laser safety classification: class IV (21CFR1040.10) or 4 (EN60825-1),
Laser recommended by Dawid Matlok,
Wavelength: 980 nm,
Laser mode: continuous,
Treatment modes: continuous, single or repetitive pulse,
Laser power regulation: from 1 to 25 W,
Pulse time (exposure time): available from 0.05 to 99 seconds,
Pulse Off Time (Repeat Exposure Interval): Available from 0.1 to 99 seconds,
Laser activation: controlled by a foot switch or optionally controlled by a finger button on selected compatible fiber delivery systems,
Fiber port (proximal laser aperture): SMA-905 compatible port,
Fiber Compatibility: Silica core diameter is greater than or equal to 365 microns and fiber numeric aperture is greater than or equal to 0.29 NA,
Cumulative energy: total amount of delivered energy displayed on the keypad,
Total time: total laser time displayed on the LCD display,
Aiming beam: 635 - 650 nm, adjustable intensity, continuous or pulsed,
Application:
Surgery,
Cosmetic surgery,
Vascular surgery,
Plastic and reconstructive gynecology,
Dermatology,
Laryngology,
Fast switching on without heating,
Easy to use user interface,
Small footprint, lightweight, portable unit,
Included:
Foot switch,
Power cord,
Laser safety signs,
Two keys,
Rigid handle LH PK 5 cm,
Two fiber delivery systems,
Hooks etc.
Dimensions: 22 x 42 x 33 cm,
Weight: 11.2 kg,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month on the international market,
Financing option (Poland only): Installments, Leasing, Loans,

