



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion pump (syringe) GRASEBY 3100 (Reconditioned)

Used, good technical condition,
English production,
Power supply: 230 V and via a built-in battery,
Frequency: 50/60 Hz,
Power consumption: 15 VA,
Reliable pump for infusion of painkillers,
Remarkable ease of use,
Durable construction,
Choice of infusion units,
Unique intermittent infusion mode
Preselection of the infusion volume or time,
Programmable in mass units,
A wide range of syringe types and sizes,
Event registration,
Communication via the RS232 port,
Foil control panel and liquid crystal display,
Acoustic alarms,
Syringe selection: 5, 10, 20, 30, 50/60 ml,
Infusion rate range (for a 50 ml syringe): 0.1 - 199.9 ml / h,
Included:
Hanging handle,
Power cord,
Dimensions: 32.5 x 22.5 x 11 cm,
Weight: 3 kg,
User manual in English (PDF),
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month on the international market,
Financing option (Poland only): Installments, Leasing, Loans,
Prześlij opinię
Historia
Zapisane
Społeczność

