



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CLV-S30 light source (Reconditioned)

Used, very good technical condition,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Adapted to rigid endoscopes,
Xenon and halogen light,
Source of light:
300 W xenon bulb - work lighting,
150 W halogen bulb - emergency lighting (possibility to work with halogen lighting),
5-step light intensity adjustment,
Color temperature: 6000K,
Dimensions: 36 x 40 x 13.5 cm,
Weight: 12.3 kg,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Installments, Leasing, Loan,

