



GE DINAMAP PRO 300V2 patient monitor (Reconditioned)

Used, good technical condition,

U.S.A. production,

Power supply: 230 V and via a 12V / 2.3Ah battery (new battery included),

Frequency: 50/60 Hz,

Power consumption: 30 VA,

The DINAMAP PRO 300V2 device by GE is used for non-invasive and automatic measurement of systolic and diastolic blood pressure, heart rate and oxygen saturation (SpO2) in adult and pediatric patients,

The device can be used in a wide variety of healthcare facilities. These include hospital departments as well as other types of medical facilities, such as doctor's offices, clinics and long-term care facilities,

Large readable liquid crystal display and smaller digital displays,

Simple and intuitive operation,

Built-in printer with a paper width: 60 mm,

Cuff pressure range:

0 - 290 mmHg (adult),

0 - 140 mmHg (newborn),

Pulse rate range: 30-200 beats per minute

SpO2 measurement range: 1 - 100%,

Included:

Air Tube Cuff,

BIOGENESIS SpO2 silicone sensor,

Printer paper roll,

Power cord,

Dimensions: 25 x 17.5 x 25 cm,

Weight: 3.5 kg,

User manual in English (PDF),

Cat 02,

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Installments, Leasing, Loan,

